**Biosigma s.r.l.**

**a Dominique Dutscher Company**

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |

Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | http://www.biosigma.com |

**M521D**

Rev. 09/03/2016

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**USE**

Open the peel pouch as indicated. Take the swab by the cap. Collect the specimen with swab's bud. Open the tube and place the swab into the gel until the cap closes firmly the tube.

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**FEATURES**

STUAR Transport Medium prepared according to the Stuart formula, is a semisolid, highly reducing, non-nutritive medium used for the transport of specimens of clinical source to laboratory. Stuart's preparation was designed to preserve Neisseriae Gonorrhoeae and Trichomonas Vaginalis without permitting their proliferation. The medium has been shown to be suitable for transporting Haemophilus Influentiae, streptococcus Pyogenes, Corynebacterium Diptheriae, Enterobacteria and Bacteria of the respiratory tract as well. Bacteria can survive up to 72 hours (resistant species) but it is better to make analysis as soon as possible

**FORMULA (grammi per litro)**

SodiumThyoglicollate.....	0,90
Sodium glycerophosphate.....	10
Calcium Chloride.....	0,10
Agar bios.....	4
pH = 7.2 ± 0.2	

**FORMULA (grams per litre)**

SodiumThyoglicollate.....	0,90
Sodium glycerophosphate.....	10
Calcium Chloride.....	0,10
Agar bios.....	4
pH = 7.2 ± 0.2	

**BIBLIOGRAPHY**

Stuart, R.D. (1954), The problem of transport of specimens for culture of Gonococchi. Can. J. Pubbl. Health, 45,73-83.

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AMIES TRANSPORT MEDIUM *without CHARCOAL (CHIARO)*

Terreno per raccolta e trasporto materiali clinici.
Medium for collection and transport of clinical samples.

BSM003

Table with 3 columns: PRODOTTO / PRODUCT, CODICE / CODE, CONFEZIONE / PACKAGING. Rows describe different tamponi and their packaging sizes.

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FEATURES

AMIES TRANSPORT MEDIUM w/o CHARCOAL is a semi-solid medium, non nutritional, utilized for the transport of the specimens to submit to bacteriology analysis. The formula is a modification of Stuart's medium containing sodium glycerophosphate utilized as a source of energy by coliforms. The thioglycollate sodium lowers the potential redox of the medium, allowing a better preservation of the anaerobic bacteria. The presence of the buffers prevents the overgrowth of the banal contaminants. The medium allows good conservation of beta-hemolytic streptococci, intestinal pathogenic and S. aureus for over three days, of N. gonorrhoeae for 24 hours, of B. pertussis and Fusobacterium for 48 to 72 hours.

FORMULA (grammi per litro)

Table listing ingredients and their amounts for the formula in grams per liter.

FORMULA (grams per litre)

Table listing ingredients and their amounts for the formula in grams per liter.

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Table showing growth characteristics for Haemophilus influenzae ATCC 10211 on Chocolate Agar.

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
**FEATURES**

AMIES TRANSPORT MEDIUM is a semi-solid medium, non nutritional, utilized for the transport of the specimens to submit to bacteriology analysis. The formula is a modification of Stuart's medium containing sodium glycerophosphate utilized as a source of energy by coliforms. The vegetable carbon present in the medium improves the survival of *Neisseria gonorrhoeae* and of the intestinal pathogenic bacteria neutralizing the toxic metabolic products. The thioglycollate sodium lowers the potential redox of the medium, allowing a better preservation of the anaerobic bacteria. The presence of the buffers prevents the overgrowth of the banal contaminants. The medium allows good conservation of beta-hemolytic streptococci, intestinal pathogenic and *S. aureus* for over 3 days, of *N. gonorrhoeae* for 24 hours, of *B. pertussis* and *Fusobacterium* for 48 - 72 hours.

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