

Code No.12-0672

SPECIFICATION

Lysozyme, from Egg White for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, powder or granules
Solubility in water	to pass test
Specific absorbance E1%1cm (λ max)	under investigation
Loss on drying at 105°C	max.8.0%
pH of a 15g/L solution at 25°C	3.0~5.0
Chloride (Cl)	max.4.2%
Potency (calculated on the dried basis)	min.0.8mg Lysozyme/mg