

Code No.14-0401

SPECIFICATION

Nitric Acid (1.38)

for Analysis of Poisonous Metals

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Zinc (Zn)	max.5ppb
Cadmium (Cd)	max.0.1ppb
Mercury (Hg)	max.0.05ppb
Tin (Sn)	max.10ppb
Lead (Pb)	max.0.5ppb
Arsenic (As)	max.1ppb
Antimony (Sb)	max.10ppb
Vanadium (V)	max.10ppb
Selenium (Se)	max.10ppb
Chromium (Cr)	max.1ppb
Manganese (Mn)	max.1ppb
Assay	60.0~62.0%