

Code No.08-0107

SPECIFICATION

Hydrochloric Acid

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Assay	35.0~37.0%
Sulfate (SO ₄)	max.3ppm
Free chlorine	to pass test (Cl:max.about 0.2ppm)
Heavy metals (as Pb)	max.2ppm
Arsenic (As)	max.0.1ppm
Iron (Fe)	max.1ppm
Substances reducing iodine	to pass test (as SO ₃ :max.about 4ppm)