

SPECIFICATION

Sodium Azide

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ nearly white, crystalline powder or
Assay (NaN ₃)	min.98.0%(mass/mass)
Solubility in water	to pass test
Alkalinity (as NaOH)	max.0.2%(mass/mass)
Loss on drying at 105°C	max.0.1%(mass/mass)
Chromium (Cr)	max.2ppm(mass/mass)

JIS K 9501 : 2019(*:Additional test performed by Wako but not required by JIS)