

Code No.04-0721

SPECIFICATION

Dimethyl Sulfoxide

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay [(CH ₃) ₂ SO](GC)	min.99.0%(area/area)
Density (20°C)	1.099~1.103g/mL
Freezing point	min.18°C
Refractive index 20°C	1.477~1.480
Water	max.0.2%(mass/mass)

JIS K 9702 : 2025(*:Additional test performed by Wako but not required by JIS)