

Code No.02-0042

SPECIFICATION

Barium Sulfate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White powder
Precipitation volume	min.2.2mL/g
Loss on ignition at 700°C	max.2.0%
Neutrality of aqueous slurry	to pass test
Phosphate	to pass test
Sulfide	to pass test
Soluble barium salts	to pass test
Heavy metals (as Pb)	max.0.001%
Arsenic (As)	max.1ppm
Iron (Fe)	max.0.005%