

Code No.13-1829

SPECIFICATION

Methyl Parahydroxybenzoate

for the Japanese Pharmacopoeia General Tests

REQUIREMENT	SPECIFICATION
Description	Colorless crystals or white crystalline powder
Identification	to pass test
Melting point	125~128°C
Purity	-
(1) Clarity and color of solution	to pass test
(2) Acidity	to pass test
(3) Related substances	to pass test
Residue on ignition	not more than 0.1%
Assay	98.0~102.0%
* Loss on drying over silica gel in vacuum	max.0.5%

*:Additional test performed by Wako