

Code No.02-1873

SPECIFICATION

Butyl Parahydroxybenzoate

for the Japanese Pharmacopoeia General Tests

REQUIREMENT	SPECIFICATION
Description	Colorless crystals or white crystalline powder
Identification	to pass test
Melting point	68~71℃
Purity	-
(1) Clarity and color of solution	to pass test
(2) Acidity	to pass test
(3) Related substances	to pass test
Residue on ignition	not more than 0.1%
Assay	98.0~102.0%
* Loss on drying over silica gel in vacuum	under investigation

*:Additional test performed by Wako