

SAFETY DATA SHEET

According to JIS Z 7253:2019
Revision date 16-Dec-2024
Revision Number 1.07

Section 1: PRODUCT AND COMPANY IDENTIFICATION

Product Name	Orexin A ELISA Kit Wako
Product Code	293-79801

Supplier FUJIFILM Wako Pure Chemical Corporation
1-2 Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan
Phone: +81-6-6203-3741
Fax: +81-6-6203-2029

Emergency telephone number +81-6-6203-3741 / +81-3-3270-8571

Recommended uses For research use only

Restrictions on use Seek expert judgment when using for purposes other than those recommended.

Section 2: HAZARDS IDENTIFICATION

GHS classification

Classification of the substance or mixture

Corrosive to metals

Category 1

Acute toxicity - Inhalation (Dusts/Mists)

Category 3

Skin corrosion/irritation

Category 1

Serious eye damage/eye irritation

Category 1

Skin sensitization

Category 1

Specific target organ toxicity (single exposure)

Category 1

Category 1 central nervous system, blood system, respiratory system

Specific target organ toxicity (repeated exposure)

Category 1

Category 1 central nervous system, respiratory system

Acute aquatic toxicity

Category 3

Chronic aquatic toxicity

Category 2

Pictograms



Signal word

Danger

Hazard statements

H290 - May be corrosive to metals

H314 - Causes severe skin burns and eye damage

H318 - Causes serious eye damage

H331 - Toxic if inhaled

H317 - May cause an allergic skin reaction

H402 - Harmful to aquatic life

H411 - Toxic to aquatic life with long lasting effects

H370 - Causes damage to the following organs: central nervous system, blood system, respiratory system

H372 - Causes damage to the following organs through prolonged or repeated exposure: central nervous system, respiratory system

Precautionary statements-(Prevention)

- Use only outdoors or in a well-ventilated area
- Do not breathe dust/fume/gas/mist/vapors/spray
- Wash face, hands and any exposed skin thoroughly after handling
- Wear protective gloves/protective clothing/eye protection/face protection
- Contaminated work clothing should not be allowed out of the workplace
- Do not eat, drink or smoke when using this product
- Avoid release to the environment
- Keep only in original container

Precautionary statements-(Response)

- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
- Immediately call a POISON CENTER or doctor/physician
- IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower
- Wash contaminated clothing before reuse
- If skin irritation or rash occurs: Get medical advice/attention
- IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing
- IF SWALLOWED: Rinse mouth. Do NOT induce vomiting
- Collect spillage
- Absorb spillage to prevent material damage

Precautionary statements-(Storage)

- Store in a well-ventilated place. Keep container tightly closed
- Store locked up
- Store in corrosive resistant/ container with a resistant inner liner

Precautionary statements-(Disposal)

- Dispose of contents/container to an approved waste disposal plant

Others

Other hazards Not available

Section 3: COMPOSITION/INFORMATION ON INGREDIENTS

Single Substance or Mixture Kit (Set of mixtures)

Chemical Name	Weight-%	Molecular weight	ENCS	ISHL No.	CAS RN
Antibody-coated 96-well Plate	-	N/A	N/A	N/A	N/A-29-7981
Orexin A Standard	-	N/A	N/A	N/A	N/A-29-7982
Anti Orexin A, Biotin-conjugated	-	N/A	N/A	N/A	N/A-29-7983
Peroxidase-conjugated Streptavidin Solution	-	N/A	N/A	N/A	N/A-29-7984
TMB Solution	-	N/A	N/A	N/A	N/A-29-7985
Buffer	-	N/A	N/A	N/A	N/A-29-7986
Stop Solution	-	N/A	N/A	N/A	N/A-29-7987
Wash Solution(20×)	-	N/A	N/A	N/A	N/A-29-7988
Plate Seal	-	N/A	N/A	N/A	N/A-29-7989

Note on ISHL No.: * in the table means announced chemical substances.

Substances Remarks:

This Product includes the following componets. Sulfuric Acid <10 %, Polyethylene glycol p-octylphenyl ether <5.0 %, 1,2-Benzisothiazol-3(2H)-one <0.50 %, 1,2-Propanediol <2.0 %, Sodium Chloride 60 - 70 %

Section 4: FIRST AID MEASURES

Inhalation

Remove to fresh air. If symptoms persist, call a physician.

Skin contact

Wash off immediately with soap and plenty of water. If symptoms persist, call a physician.

Eye contact

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediate medical attention is required.

Ingestion

Rinse mouth. Never give anything by mouth to an unconscious person. Call a physician or poison control center immediately. Do not induce vomiting without medical advice.

Protection of first-aiders

Use personal protective equipment as required.

Section 5: FIRE FIGHTING MEASURES**Suitable extinguishing media**

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment

Unsuitable extinguishing media

No information available

Specific hazards arising from the chemical product

Thermal decomposition can lead to release of irritating and toxic gases and vapors.

Special extinguishing method

No information available

Special protective actions for fire-fighters

Use personal protective equipment as required. Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear.

Section 6: ACCIDENTAL RELEASE MEASURES**Personal precautions, protective equipment and emergency procedures**

For indoor, provide adequate ventilation process until the end of working. Deny unnecessary entry other than the people involved by, for example, using a rope. While working, wear appropriate protective equipments to avoid adhering it on skin, or inhaling the gas. Work from windward, and retract the people downwind.

Environmental precautions

To be careful not discharged to the environment without being properly handled waste water contaminated.

Methods and materials for contaminant and methods and materials for cleaning up

Absorb dry sand, earth, sawdust and the waste. Collect empty container that can be sealed.

Recovery, neutralization

No information available

Secondary disaster prevention measures

Clean contaminated objects and areas thoroughly observing environmental regulations.

Section 7: HANDLING AND STORAGE**Handling****Technical measures**

Avoid contact with alkaline substances. Avoid contact with metal. Use with local exhaust ventilation.

Precautions

Do not rough handling containers, such as upsetting, falling, giving a shock, and dragging. Prevent leakage, overflow, and scattering. Not to generate steam and dust in vain. Seal the container after use. After handling, wash hands and face, and then gargle. In places other than those specified, should not be smoking or eating and drinking. Should not be brought contaminated protective equipment and gloves to rest stops. Deny unnecessary entry of non-emergency personnel to the handling area.

Safety handling precautions

Use personal protective equipment as required. Avoid contact with skin, eyes or clothing.

Storage**Safe storage conditions****Storage conditions**

Store away from sunlight in a cool (2-10 °C) well-ventilated dry place.

Safe packaging material

No information available

Incompatible substances

alkaline substances, Metals

Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

Engineering controls

In case of indoor workplace, seal the source or use a local exhaust system. Provide the safety shower facility, and hand- and eye-wash facility. And display their position clearly.

Exposure limits

Chemical Name	JSOH (Japan)	ISHL (Japan)	ACGIH
Sulfuric Acid 7664-93-9	Ceiling: 1 mg/m ³	N/A	TWA 0.2mg/m ³

Personal protective equipment

Respiratory protection	Gas mask for acidic gas (JIS T 8152)
Hand protection	chemical protective gloves (JIS T 8116)
Eye protection	protective eyeglasses or chemical safety goggles (JIS T 8147)
Skin and body protection	Long-sleeved work clothes

General hygiene considerations

Handle in accordance with good industrial hygiene and safety practice.

If this product is classified as "Chemical Substances Hazardous to Skin, etc.", use appropriate protective equipment to them.

Section 9: PHYSICAL AND CHEMICAL PROPERTIES

Form

Appearance	Kit (Set of mixtures)
Odor	no data available
Melting point/freezing point	no data available
Boiling point, initial boiling point and boiling range	no data available
Flammability	no data available
Evaporation rate:	no data available
Flammability (solid, gas):	no data available
Upper/lower flammability or explosive limits	
Upper:	no data available
Lower:	no data available
Flash point	no data available
Auto-ignition temperature:	no data available
Decomposition temperature:	no data available
pH	no data available
Viscosity (coefficient of viscosity)	no data available
Dynamic viscosity	no data available
Solubilities	no data available
n-Octanol/water partition coefficient:(log Pow)	no data available
Vapour pressure	no data available
Specific Gravity / Relative density	no data available
Vapour density	no data available
Particle characteristics	no data available

Section 10: STABILITY AND REACTIVITY

Stability

Reactivity	no data available
Chemical stability	Stable under recommended storage conditions.
Hazardous reactions	
	Corrodes metals to generate hydrogen gas.
Conditions to avoid	

Extremes of temperature and direct sunlight

Incompatible materials

alkaline substances, Metals

Hazardous decomposition products

Carbon monoxide (CO), Carbon dioxide (CO₂), Nitrogen oxides (NO_x), Sulfur oxides (SO_x)

Section 11: TOXICOLOGICAL INFORMATION

*NITE: National Institute of Technology and Evaluation (JAPAN)

https://www.chem-info.nite.go.jp/en/chem/chrip/chrip_search/srhInput

Acute toxicity

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Sulfuric Acid	2140 mg/kg (Rat)	N/A	0.375 mg/L (Rat) 4 h
1,2-Propanediol	22000 mg/kg (Rat)	20800 mg/kg (Rabbit)	N/A
Polyethylene glycol p-octylphenyl ether	1800 mg/kg (Rat)	N/A	N/A
1,2-Benzothiazolin-3-one	670 mg/kg (Rat)	> 2000 mg/kg bw (Rat)	0.25 mg/L (Rat)

Chemical Name	Acute toxicity -oral- source information	Acute toxicity -dermal- source information	Acute toxicity -inhalation gas- source information
Sulfuric Acid	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.
1,2-Propanediol	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.

Chemical Name	Acute toxicity -inhalation vapor- source information	Acute toxicity -inhalation dust- source information	Acute toxicity -inhalation mist- source information
Sulfuric Acid	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.
1,2-Propanediol	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.

Skin irritation/corrosion

Chemical Name	Skin corrosion/irritation source information
Sulfuric Acid	Based on the NITE GHS classification results.
1,2-Propanediol	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.

Serious eye damage/ irritation

Chemical Name	Serious eye damage/irritation source information
Sulfuric Acid	Based on the NITE GHS classification results.
1,2-Propanediol	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.

Respiratory or skin sensitization

Chemical Name	Respiratory or Skin sensitization source information
Sulfuric Acid	Based on the NITE GHS classification results.
1,2-Propanediol	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.

Reproductive cell mutagenicity

Chemical Name	germ cell mutagenicity source information
Sulfuric Acid	Based on the NITE GHS classification results.

1,2-Propanediol	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.

Carcinogenicity

Chemical Name	Carcinogenicity source information
Sulfuric Acid	Based on the NITE GHS classification results.
1,2-Propanediol	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.

Chemical Name	NTP	IARC	ACGIH	JSOH
Sulfuric Acid 7664-93-9	-	Group 1	A2	-

Reproductive toxicity

Chemical Name	Reproductive toxicity source information
Sulfuric Acid	Based on the NITE GHS classification results.
1,2-Propanediol	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.

STOT-single exposure

Chemical Name	STOT -single exposure- source information
Sulfuric Acid	Based on the NITE GHS classification results.
1,2-Propanediol	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.

STOT-repeated exposure

Chemical Name	STOT -repeated exposure- source information
Sulfuric Acid	Based on the NITE GHS classification results.
1,2-Propanediol	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.

Aspiration hazard

Chemical Name	Aspiration Hazard source information
Sulfuric Acid	Based on the NITE GHS classification results.
1,2-Propanediol	Based on the NITE GHS classification results.
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results.
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results.

Section 12: ECOLOGICAL INFORMATION

*NITE: National Institute of Technology and Evaluation (JAPAN)
https://www.chem-info.nite.go.jp/en/chem/chrip/chrip_search/srhInput

Ecotoxicity

Chemical Name	Algae/aquatic plants	Fish	Crustacea
Sulfuric Acid	N/A	LC50: <i>Lepomis macrochirus</i> 16 - 28 mg/L 96 h	LC50: <i>Daphnia magna</i> 29 mg/L 24 h
1,2-Propanediol	EC50 : <i>Pseudokirchneriella subcapitata</i> > 1000 mg/L 72 h	LC50 : <i>Oryzias latipes</i> > 100 mg/L 96 h	EC50 : <i>Daphnia magna</i> > 1000 mg/L 48 h
Polyethylene glycol p-octylphenyl ether	N/A	LC50 : <i>Lepomis macrochirus</i> 3 mg/L 96 h	N/A
1,2-Benzothiazolin-3-one	ErC50 : <i>Pseudokirchneriella subcapitata</i> 0.22 mg/L 72 h	N/A	N/A

Other data

Chemical Name	Short-term (acute) hazardous to the	Long-term (chronic) hazardous to the
---------------	-------------------------------------	--------------------------------------

	aquatic environment source information	aquatic environment source information
Sulfuric Acid	Based on the NITE GHS classification results	Based on the NITE GHS classification results
1,2-Propanediol	Based on the NITE GHS classification results	Based on the NITE GHS classification results
Polyethylene glycol p-octylphenyl ether	Based on the NITE GHS classification results	Based on the NITE GHS classification results
1,2-Benzothiazolin-3-one	Based on the NITE GHS classification results	Based on the NITE GHS classification results

Persistence and degradability	No information available
Bioaccumulative potential	No information available
Mobility in soil	No information available
Hazard to the ozone layer	No information available

Section 13: DISPOSAL CONSIDERATIONS

Waste from residues

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Contaminated container and contaminated packaging

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Section 14: TRANSPORT INFORMATION

ADR/RID

UN number	UN2796
Proper shipping name:	Sulphuric acid
UN classification	8
Subsidiary hazard class	
Packing group	II
Marine pollutant	Yes

IMDG

UN number	UN2796
Proper shipping name:	Sulphuric acid
UN classification	8
Subsidiary hazard class	
Packing group	II
Marine pollutant (Sea)	Yes
Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	No information available

IATA

UN number	UN2796
Proper shipping name:	Sulphuric acid
UN classification	8
Subsidiary hazard class	
Packing group	II
Environmentally Hazardous Substance	Yes

Section 15: REGULATORY INFORMATION

Japanese regulations

Fire Service Act	Not applicable
Poisonous and Deleterious Substances Control Law	Not applicable
Industrial Safety and Health Act	Harmful Substances Whose Names Are to be Indicated on the Label (Law Art.57)

Industrial Safety and Health Act (2025~)	Notifiable Substances (Law Art.57-2) Group 3 Specified Chemical Substance, (Ordinance on Prevention of Hazards Due to Specified Chemical Substances Art.2 Para.1, Item 6) Chemical Substances Hazardous to Skin, etc.(Regulations Article 594-2 Paragraph 1) 【2025.4.1~】Harmful Substances Whose Names Are to be Indicated on the Label (Law Art.57) 【2025.4.1~】Notifiable Substances (Law Art.57-2)
Industrial Safety and Health Act (2026~)	【2026.4.1~】Harmful Substances Whose Names Are to be Indicated on the Label (Law Art.57) 【2026.4.1~】Notifiable Substances (Law Art.57-2)
Regulations for the carriage and storage of dangerous goods in ship	Corrosive Substances (Ordinance Art.3, Ministry of Transportation Ordinance Regarding Transport by Ship and Storage, Attached Table 1)
Civil Aeronautics Law	Corrosive Substances (Ordinance Art.194, MITL Notification for Air Transportation of Explosives etc., Attached Table 1)
Marine Pollution Prevention Law	Enforcement ordinance Appendix No. 1 Noxious liquid substance Category Y
Pollutant Release and Transfer Register Law (2023.4.1-)	Class 1
Class 1 - No.	408
Water Pollution Control Act	Specified substances(Law Art.2 Para.4, Enforcement Order Art.3-3)
Export Trade Control Order	Not applicable
Air Pollution Control Law	Specified Substances

Pollution Release and Transfer Registry (~2023.3.31)

Class	Chemical Name in Regulation	(Metal Name)	Control number	Content Rate
Class 1	Poly(oxyethylene) octylphenyl ethers	-	408	<5.0

Industrial Safety and Health Law

Law Name	Chemical Name in Regulation	Weight %	
Notifiable Substances (Law Art.57-2)	Sulfuric acid	<10	Existing Law
Notifiable Substances (Law Art.57-2)	(±)-Propylene glycol	<2.0	2025/4/1
Notifiable Substances (Law Art.57-2)	1,2-Benzisothiazol-3(2H)-one	<0.50	2025/4/1
Notifiable Substances (Law Art.57-2)	Poly(oxyethylene) p-(1,1,3,3-tetramethylbutyl)phenyl ether	<5.0	2026/4/1

Section 16: OTHER INFORMATION**Key literature references and sources for data etc.**

NITE: National Institute of Technology and Evaluation (JAPAN)
https://www.chem-info.nite.go.jp/en/chem/chrip/chrip_search/srhInput
 IATA dangerous Goods Regulations
 RTECS:Registry of Toxic Effects of Chemical Substances
 Japan Industrial Safety and Health Association GHS Model SDS
 Dictionary of Synthetic Organic Chemistry, SSOCJ, Koudansha Scientific Co.Ltd.
 Chemical Dictionary, Kyouritsu Publishing Co., Ltd.
 etc

Record of SDS revisions

The following contents were revised. Hazards identification. Composition/information on ingredients. Accidental release measures. Handling and storage. Stability and reactivity. Toxicological information. Ecological information. Transport information. Regulatory information.

Disclaimer

This SDS is according to JIS Z 7253: 2019. The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

GHS Classification is according to JIS Z 7252:2019. *JIS: Japanese Industrial Standards

End of Safety Data Sheet