

Code No.01-0449

SPECIFICATION

Antimony(III) Chloride

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ slightly yellow, crystals ~ crystalline powder
Assay	min.98.0%
Solubility in dil. hydrochloric acid	to pass test
Sulfur compounds (as SO ₄)	max.0.01%
Sodium (Na)	max.0.005%
Copper (Cu)	max.5ppm
Lead (Pb)	max.0.003%
Arsenic (As)	max.0.002%
Iron (Fe)	max.0.002%

JIS K 8400-1994(*:Additional test performed by Wako but not required by JIS)