

Code No.01-1860

## SPECIFICATION

### Acequinocyl-hydroxy Standard

for Pesticide Residue Analysis

| REQUIREMENT                              | SPECIFICATION             |
|------------------------------------------|---------------------------|
| Appearance                               | Yellow crystalline powder |
| Solubility in dichloromethane            | to pass test              |
| Melting point                            | 94~98°C                   |
| Loss on drying over silica gel in vacuum | under investigation       |
| Assay (HPLC)                             | min.99.0%                 |
| Assay (qNMR)                             | min.99.0%                 |